

Sanidad Facial Toner



MODEL: RP 8411

- ▶ Tender 1.7Hz AcuPulse-TENS; 1400Gauss Magnetic Energy.
- ▶ Wireless Sponge-Pad Output; Fully Moisturized, good for a day's use.
- ▶ Easy DIY Acupuncture Therapy Manual; Safety Water-Proof Design

磁波
、
針灸
、
美容

~ Your Pocket Magnetic Acu-Pulser ~

MEDICAL-APPLICATION PROPERTIES:

- ⊕ A fixed 1.7Hz low frequency TENS treatment.
- ⊕ Non-invasive, drug-free method of pain relief.
- ⊕ The body's natural endorphins are triggered as analgesics after 20 minutes of application. (drink water after treatment)
- ⊕ Treating both acute and chronic pains.
- ⊕ Facilitating blood circulation and muscle exercises (when muscle motor points are triggered).
- ⊕ The web sponge electrodes provide tender touch in TENS sensation.



SPECIFICATIONS:

- Main device diameter/height: 5.4cm/2.3cm (2.1inch/0.90inch)
- 1.7Hz Low Frequency Acu-pulse; 1400G Magnet x 2
- DC 6V, 4 bottom Alkali batteries; Power Indicator: Red Light
- 2 sponge electrode pads (PVA-3315); 1 protection cover cap.
- 1 Magnet Testing tube; 1 spare pair of sponge pads, 1 screwdriver to manage battery chamber
- User manual, Gift box and 1 carrying case and leather bag.



ADDITIONAL ACCESSORIES:

- ⊕ New Patented Wire kit with 2pcs 25m/m round self-adhesive electrodes.
- ⊕ Ear-Clam inductors



WARNINGS:

- Do not apply the device under the following conditions:
(Please refer to user manual for EPI info before use)
- ⊗ Infants or children ⊗ Skin disease ⊗ Heart disease
 - ⊗ Pregnancy ⊗ Carrying cardiac pace makers, etc.



CE-M.D.D./ISO 9002 / EN 46002
 R.O.C. Medical Device
 No. 000043: Electronic Acupuncture,

Rolor Electronics Corp.

No. 274, Sec. 5, Nanking E. Rd., Taipei, Taiwan
 TEL: 886-2-27682174-6
 FAX: 886-2-27605922-3
 E-Mail: webmaster@rolor.com

