

Sanidad Facial Toner / Mini-TENS



MODEL: RP 8411N


- ▶ Tender 1~15Hz AcuPulse
Low-frequency; 1400Gauss Magnetic Energy.
- ▶ Wireless Magnetized Sponge-Pad Output
- ▶ 2 lead-wire outputs good for different types of electrodes

磁波
、
針灸
、
美容

~ Portable Magnetic Acu-Pulser ~

MEDICAL-APPLICATION PROPERTIES:

- ⊕ Adjustable 1~15Hz low-frequency acu-pulsing treatment.
- ⊕ Non-invasive, drug-free method of pain relief.
- ⊕ The body's natural endorphins are triggered as analgesics after 20 minutes of application. (drink water after treatment)
- ⊕ Treating both acute and chronic pains.
- ⊕ Facilitating blood circulation and muscle exercises (when motor nerves/muscle motor points are triggered).
- ⊕ The web sponge electrodes provide tender touch treatment with magnetic energy.



SPECIFICATIONS:

- Main device Dimension/Weight: 58x58x31m/m, 70g;
- DC 6V, 4 button Alkali batteries;
- 1~15Hz Low Frequency, 1400Gauss Magnetic Energy;
- Power Indicator: LED Green Light;
- 2 lead wires and 4 pieces of Self-adhesive Electrode pads (RUB-4040),
- 2 sponge electrode pads (PVA-2800) and 1 magnet wave testing tube;
- 1 protection cover cap, 1 spare pair of sponge pads;
- User manual, carrying leather bag, Gift-wrap box.



ADDITIONAL ACCESSORIES:

- † Ear-Clam Inductors.
- † Eye-Nose Inductor.



WARNINGS:

Do not apply the device under the following conditions: (Please refer to user manual for EPI info before use)

- ⊗ Infants or children ⊗ Skin disease ⊗ Heart disease ⊗ Pregnancy ⊗ Carrying cardiac pace makers, etc.



CE-M.D.D./ISO 9002 / EN 46002
R.O.C. Medical Device
No. 000043: Electronic Acupuncture,



Rolor Electronics Corp.

No. 274, Sec. 5, Nanking E. Rd., Taipei, Taiwan
TEL: 886-2-27682174-5
FAX: 886-2-27605922-3
E-Mail: webmaster@rolor.com
URL: <http://www.rolor.com>

